

# SPARTA DSP

## *DSP Plotting Guide: Operation Guide*

Matthew Duong (3223 Affiliate)


Created: December 13th, 2024

Last Updated: December 13th, 2024


---

### Open the Plotting Template:

1. Open the excel file called 'DSP TEMPLATE.xlsx' in the middle right side of the desktop.




2. Look at the bottom of the file and keep going to the right until you see the tabs called 'Dat0', 'Dat1', 'Dat2', ... 'DatX'. Start with 'Dat0'.




- Each of these tabs handle one DSP's primary data. So if you use five of them for example, the final plot will have five different lines on it.

### Selecting a DSP file:


- Click on 'File' at the top left corner




5. Click on 'Open'




6. Click on 'Browse'



7. Change the file type filter to 'All Files (\*.\*)' by clicking on the drop down menu



## 8. Locate and open the DSP .idf file that you want to plot



## 9. Text Import Wizard will open, click 'Next'

Text Import Wizard - Step 1 of 3

The Text Wizard has determined that your data is Delimited.

If this is correct, choose Next, or choose the data type that best describes your data.

Original data type

Choose the file type that best describes your data:

☒ Delimited - Characters such as commas or tabs separate each field.

☐ Fixed width - Fields are aligned in columns with spaces between each field.

Start import at row: 1 File origin: Windows (ANSI)

☐ My data has headers.


Preview of file C:\Users\sparta\Desktop\JSC1a\_1pct\_0084\_230620EISScan1\_P28659.idf.

```
1 |
2 | Version=12
3 | QR=QR,2,230620123053000P28659,,84,45097.521447,45097.522928,user1,100,,P28659
4 | Peaks=false
5 | Corrosion=false
6 | Levels=false
7 | ECN=false
8 | Solarcell=false
```


< >

Cancel < Back **Next >** Finish

10. Change the delimiter to check ONLY 'Space' and click next




11. On the last page, just click finish




12. This will open up a new excel window holding the DSP data.

13. It will also open up this window about converting, just close it and ignore it



## Copying DSP data to the plotter:

14. Scroll down until you see a line that is labelled 'primary\_data'. It's usually around line 140



140	primary_data			
141	3			
142	26			
143	-1.90E+03	-6.75E+04	1.00E+05	
144	-5.77E+02	-1.06E+05	6.31E+04	
145	2.38E+02	-1.65E+05	3.98E+04	
146	9.24E+02	-2.59E+05	2.51E+04	
147	1.13E+03	-4.09E+05	1.58E+04	
148	5.85E+02	-6.48E+05	1.00E+04	
149	1.62E+04	-1.02E+06	6.31E+03	
150	1.92E+04	-1.60E+06	3.98E+03	
151	-1.93E+02	-2.57E+06	2.51E+03	
152	9.71E+03	-4.06E+06	1.58E+03	
153	-4.99E+04	-6.41E+06	1.00E+03	
154	1.08E+06	-1.01E+07	6.31E+02	
155	1.16E+06	-1.60E+07	3.98E+02	
156	1.29E+06	-2.52E+07	2.51E+02	
157	1.47E+06	-3.97E+07	1.58E+02	
158	7.51E+05	-6.21E+07	1.00E+02	
159	1.42E+07	-9.30E+07	6.31E+01	
160	5.04E+06	-1.58E+08	3.98E+01	
161	-6.60E+06	-2.55E+08	2.51E+01	
162	-2.18E+07	-3.76E+08	1.58E+01	
163	1.08E+07	-6.19E+08	1.00E+01	
164	-5.35E+06	-9.19E+08	6.31E+00	
165	7.05E+08	-1.41E+09	3.98E+00	
166	-1.49E+09	-2.03E+09	2.51E+00	
167	7.19E+08	-1.21E+09	1.58E+00	
168	5.74E+08	-5.01E+09	1.00E+00	
169	2			



15. This is the primary DSP data.

16. Highlight the entire block of data and copy it

The screenshot shows a spreadsheet with a table of data. The table has 4 columns: a label column, a column with values 3 and 26, and two columns of scientific notation data. A context menu is open over the data, with the 'Copy' option highlighted by a red rectangle. The menu also includes options like 'Cut', 'Paste Options', 'Smart Lookup', 'Insert...', 'Delete...', 'Clear Contents', 'Quick Analysis', 'Filter', 'Sort', 'Get Data from Table/Rar', 'New Comment', 'New Note', 'Format Cells...', 'Pick From Drop-down Li', and 'Define Name'.

	primary_data			
40				
41		3		
42		26		
43	-1.90E+03	-6.75E+04	1.00E+0	
44	-5.77E+02	-1.06E+05	6.31E+0	
45	2.38E+02	-1.65E+05	3.98E+0	
46	9.24E+02	-2.59E+05	2.51E+0	
47	1.13E+03	-4.09E+05	1.58E+0	
48	5.85E+02	-6.48E+05	1.00E+0	
49	1.62E+04	-1.02E+06	6.31E+0	
50	1.92E+04	-1.60E+06	3.98E+0	
51	-1.93E+02	-2.57E+06	2.51E+0	
52	9.71E+03	-4.06E+06	1.58E+0	
53	-4.99E+04	-6.41E+06	1.00E+0	
54	1.08E+06	-1.01E+07	6.31E+0	
55	1.16E+06	-1.60E+07	3.98E+0	
56	1.29E+06	-2.52E+07	2.51E+0	
57	1.47E+06	-3.97E+07	1.58E+0	
58	7.51E+05	-6.21E+07	1.00E+0	
59	1.42E+07	-9.30E+07	6.31E+0	
60	5.04E+06	-1.58E+08	3.98E+0	
61	-6.60E+06	-2.55E+08	2.51E+0	
62	-2.18E+07	-3.76E+08	1.58E+0	
63	1.08E+07	-6.19E+08	1.00E+0	
64	-5.35E+06	-9.19E+08	6.31E+0	
65	7.05E+08	-1.41E+09	3.98E+0	
66	-1.49E+09	-2.03E+09	2.51E+0	
67	7.19E+08	-1.21E+09	1.58E+0	
68	5.74E+08	-5.01E+09	1.00E+0	
69		2		

17. Switch back to the 'Dat0' tab on the template file

18. Click on the box under the 'Zr' labelled column and paste the data block in

Filename: Calibri 11 A A \$ % 62023\DS

B I A

Zr Zi freq w -Zi

4		0.00E+00	0.00E+00
5		0.00E+00	0.00E+00
6		0.00E+00	0.00E+00
7		0.00E+00	0.00E+00
8		0.00E+00	0.00E+00
9		0.00E+00	0.00E+00
10		0.00E+00	0.00E+00
11		0.00E+00	0.00E+00
12		0.00E+00	0.00E+00
13		0.00E+00	0.00E+00
14		0.00E+00	0.00E+00
15		0.00E+00	0.00E+00
16		0.00E+00	0.00E+00
17		0.00E+00	0.00E+00
18		0.00E+00	0.00E+00
19		0.00E+00	0.00E+00
20		0.00E+00	0.00E+00
21		0.00E+00	0.00E+00
22		0.00E+00	0.00E+00
23		0.00E+00	0.00E+00
24		0.00E+00	0.00E+00
25		0.00E+00	0.00E+00
26		0.00E+00	0.00E+00
27		0.00E+00	0.00E+00
28		0.00E+00	0.00E+00
29		0.00E+00	0.00E+00
30			

Search the menus

Cut

Copy

**Paste Options:**

Paste Special... | >

Smart Lookup

Insert Copied Cells...

Delete...

Clear Contents

Quick Analysis

Filter >

Sort >

Get Data from Table/Range...

New Comment



New Note

3	Zr	Zi	freq		w	-Zi
4	-1.90E+03	-6.75E+04	1.00E+05		6.28E+05	6.75E
5	-5.77E+02	-1.06E+05	6.31E+04		3.96E+05	1.06E
6	2.38E+02	-1.65E+05	3.98E+04		2.50E+05	1.65E
7	9.24E+02	-2.59E+05	2.51E+04		1.58E+05	2.59E
8	1.13E+03	-4.09E+05	1.58E+04		9.96E+04	4.09E
9	5.85E+02	-6.48E+05	1.00E+04		6.28E+04	6.48E
10	1.62E+04	-1.02E+06	6.31E+03		3.96E+04	1.02E
11	1.92E+04	-1.60E+06	3.98E+03		2.50E+04	1.60E
12	-1.93E+02	-2.57E+06	2.51E+03		1.58E+04	2.57E
13	9.71E+03	-4.06E+06	1.58E+03		9.96E+03	4.06E
14	-4.99E+04	-6.41E+06	1.00E+03		6.28E+03	6.41E
15	1.08E+06	-1.01E+07	6.31E+02		3.96E+03	1.01E
16	1.16E+06	-1.60E+07	3.98E+02		2.50E+03	1.60E
17	1.29E+06	-2.52E+07	2.51E+02		1.58E+03	2.52E
18	1.47E+06	-3.97E+07	1.58E+02		9.96E+02	3.97E
19	7.51E+05	-6.21E+07	1.00E+02		6.28E+02	6.21E
20	1.42E+07	-9.30E+07	6.31E+01		3.96E+02	9.30E
21	5.04E+06	-1.58E+08	3.98E+01		2.50E+02	1.58E
22	-6.60E+06	-2.55E+08	2.51E+01		1.58E+02	2.55E
23	-2.18E+07	-3.76E+08	1.58E+01		9.96E+01	3.76E
24	1.08E+07	-6.19E+08	1.00E+01		6.28E+01	6.19E
25	-5.35E+06	-9.19E+08	6.31E+00		3.96E+01	9.19E
26	7.05E+08	-1.41E+09	3.98E+00		2.50E+01	1.41E
27	-1.49E+09	-2.03E+09	2.51E+00		1.58E+01	2.03E
28	7.19E+08	-1.21E+09	1.58E+00		9.96E+00	1.21E
29	5.74E+08	-5.01E+09	1.00E+00		6.28E+00	5.01E

19. At this point, all of the automated calculations on the right of the data block should be completed




22. To view the impedance plot, look back at the bottom tabs and keep going left until you see 'Z Mod Plot (2)'



23. For every additional plot you want to add, just keep using all of the different 'Dat' tabs and the previous steps to open DSP .idf files as before. They should all appear on the 'Z Mod Plot (2)' tab.

24. When you want to save, make sure to use 'Save As' so that the template does not get overwritten.



25. There are other tabs and different types of plots on the other ones (I think phase is included), but I don't know enough about those to give you a solid answer.